

FORESITE 20-20

INTELLIGENT QUALITY CHECK PRIMER

Steven Hansen, Sai Phalke, Jose Osvaldo Alvarez, Brett Tweedie & Rachel Momanyi

09/23/2021

Our Path Forward

WHY

we exist

We deliver products and solutions to **improve the lives of people every day.**

WHERE

we are going

We aspire to be **healthcare's most trusted partner** by building upon our scale and heritage in distribution, products and solutions, while driving growth in evolving areas of healthcare through customer insights, data and analytics, and focusing our resources on what matters most.

HOW

we will grow



VALUES

what we believe



Integrity



Inclusive



Innovative



Accountable



Mission driven

BEHAVIORS

how we act

Invites curiosity

Inspires commitment

Builds partnerships

Develops self and others

What is a defect worth?



Dissatisfaction(Angry, Upset)

Reputation

Legal/Recall

Potential for patient/user harm

Future business

\$\$\$\$\$

Introduction

2021 Medical Device Recalls

Subscribe to Email Updates

f Share

🐦 Tweet

in LinkedIn

✉ Email

🖨 Print

Medical Device Recalls

2021 Medical Device Recalls

2020 Medical Device Recalls

The list below contains recalls that were issued in 2021.

2021 Medical Device Recalls

Device Name	Date
Armstrong Medical Limited Recalls AMSORB PLUS PREFILLED G-CAN 1.0L Due to Reduced Gas Flow to Patients During Anesthesia	09/23/21
Cordis Recalls Super Torque MB Angiographic Catheter with Radiopaque Marker Bands Due to Potential for Marker Bands to Move or Dislodge	09/22/21
Medtronic Recalls Pipeline Flex Embolization Devices for Risk of Delivery System Fractures During Placement, Retrieval, or Movement of Device	09/20/21
Smiths Medical Recalls NORMOFLO Irrigation Fluid Warmers and Warming Sets Due to the Possibility of Harmful Levels of Aluminum Leaching into the Fluid Path of the Warmers	09/15/21
All Ultrasound Gels and Lotions Manufactured by Eco-Med Pharmaceutical, Inc. Recalled Due to Risk of Bacteria Contamination	09/10/21
Bio-Medical Equipment Service Co. Recalls Alaris Infusion Pump Module 8100 Bezel Due to Possible Cracked or Separated Bezel Repair Posts	08/24/21
Cardinal Health Recalls Monoject Saline Flush Prefilled Syringes for Risk of Air Re-entering Syringe Leading to Air Embolism	08/23/21
Cardinal Health Recalls Argyle UVC Insertion Tray Due to Missing Instructions for Use for the Safety Scalpel N11	08/20/21

Content current as of:
09/23/2021

Regulated Product(s)
Medical Devices
Radiation-Emitting Products

Cardinal Health Recalls Monoject Saline Flush Prefilled Syringes for Risk of Air Re-entering Syringe Leading to Air Embolism



Share



Tweet



LinkedIn



Email



Print

Medical Device Recalls

[2021 Medical Device Recalls](#)

[2020 Medical Device Recalls](#)

The FDA has identified this as a Class I recall, the most serious type of recall. Use of these devices may cause serious injuries or death.

Recalled Product

- Monoject Flush Prefilled Syringes (0.9% sodium chloride)
- Product Codes, All Lots:
 - **Product Code (SKU) - Product Description:**
 - 8881570121 - 12mL Syringe, 10mL Saline Fill
 - 8881570123 - 12mL Syringe, 3mL Saline Fill
 - 8881570125 - 12mL Syringe, 5mL Saline Fill
- Distribution Dates: July 1, 2019 to July 1, 2021

Content current as of:
08/23/2021

Introduction

Problem statement:

Addressing quality issues in manufacturing during or after the production is very time consuming, expensive and adds to restrain on backorders.

Current cost of awaiting Disposition for month of August 2021 at Crystal Lake Facility: **\$ 441,517**

Average product disposition cost per year : **\$ 1,103,762**

(excluding labor, overhead costs etc)


Risk Associated: Regulatory and Compliance Risks.

*Inherent friction between quality and manufacturing.

Solution :

Leverage LOT data, maintenance data and other manufacturing inputs to predict potential quality issues so they can be addressed before production.

Manufacturing Process Overview



Data Mining



DEVICE HISTORY RECORD
CONTAINERS / CRYSTAL LAKE

Product Code: 6906 Lot #: 20406063

Start-Up: Shift: C Date: 05-12-20 End: Shift: C Date: 05-19-20

☐ N/A (I) If not applicable

	<input type="checkbox"/> N/A	<input checked="" type="checkbox"/> N/A	<input type="checkbox"/> N/A	<input type="checkbox"/> N/A	<input type="checkbox"/> N/A	<input type="checkbox"/> N/A	<input type="checkbox"/> N/A
Press # / Location	<u>70</u>					<u>70</u>	<u>70</u>
Part # with tool letter	<u>70270Y</u>	R	R	R	WP	<u>7026</u>	<u>6906</u>
Quantity Produced	<u>76449</u>						<u>756</u>

Auxiliary Equipment Number (example: heat staker, labeler, welder, rotor machine, printer etc.): 05/19/20

☐ N/A (I) If Equipment # not applicable

ID# AL-65 ID# LP-135 ID# *N/A ID# *N/A ID# *N/A ID# *N/A

Component Identification *N/A's CCoHes

☐ N/A (I) If not applicable


N/A Component	N/A Lot#	N/A Insp. By	N/A Component	N/A Lot#	N/A Insp. By
<u>Press 6610</u>	<u>1187-61</u>	<u>S. N. J. J. J.</u>			
<u>Part 15</u>	<u>1022-139</u>	<u>S. N.</u>			
<u>CL6245-139</u>	<u>*N/A</u>	<u>S. N.</u>			

Review and Approval Checklist

Checklist	Forms Presented and Approved (Check one)		Verified By:
Production Quantity Verification	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> N/A	<u>CCoHes</u>
PCRs/Start-up Checklist / Line Clearance Checklist	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> N/A	<u>CC</u>
Traceability of Components	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> N/A	<u>CC</u>
SU/PI/ER/ECR Inspection Results	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> N/A	<u>CC</u>
Reject/Rework/Sorting Forms	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> N/A	<u>CC</u>
Full Device History Record requirements have been completed and met (Print/Sign/Date)	<u>Cindy Cost</u>		<u>05/19/20</u>
NC-IL081- closed in SmartSolve?	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> N/A	<u>A Castaneda</u>
BPCS Hold closed?	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> N/A	<u>AC</u>
Non-standard shop order disposition requirements met?	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> N/A	<u>AC</u>
Full Device History Record requirements have been completed and met (Print/Sign/Date)	<u>A Castaneda</u>		<u>05-20-20</u>

FM20006690


Data Model Overview



- 10% of training data had NCs
- Class imbalance
- Small dataset
- Focused on ensemble method

----- AdaBoost low estimators, deep depth Model -----
percent of lots indicted = 1.211%
recall (NC reduction) = 9.43%
total number of predicted NCs = 7
True Positives = 5.0
precision = 62.5%

FORESITE 20-20



Impact

Process Efficiencies

- Earlier interventions on lots predicted to have a higher non-conformance rate
- Customized quality check process
- Resource efficiency
- Decreased production turnaround time.

Material Savings

- Reduce Average product cost per year by 10% (\$ **110,376**) at Crystal Lake
- Additional Savings for labor and overhead cost
- Reduction of cost of awaiting Disposition monthly at Crystal Lake: \$ **44,151**



Growth Opportunities

- Model can scale to other sites and within Cardinal
- Enhanced the **partnership** between Manufacturing and Quality

Questions

